

Cold Box

Model No: **GKS-CB-143**

Conforming to **WHO/PQS/E004**

Product type-	Short Range	Cold Life time-	More than 50 hrs
Type of Ice packs-	0.3 L	Ambient zone -	+43°C



Specifications

External dimensions (LxWxH)Approx cm	48x44x40
No of Ice packs	14
Type of Ice packs	0.3 Liters
Vaccine Storage dimensions (LxWxH) (cm)	22x22x16
Cold life at +43°C	>50 Hrs
Vaccine storage capacity	>6 Liters

Insulation Material	Polyurethane
Insulation thickness	60-75 mm
Empty Weight(kg)	10.9
Fully Loaded weight (Kg)	20
External Surface material	LLDP
Internal Lining Material	LLDP

Address

GKS Healthsol LLP
B-1/03 Palm Grove Villas, Ardee City,
Sector-52 Gurgaon Haryana,
India,122011

Factory

GKS Healthsol LLP
Plot No-362, sec-8 IMT Manesar.
Gurugram,122051 , HR,

Contact us

email:- dinesh@gkshealth.com
Mb.No:- +91-9899632558/8010743888

ISO 9001,13485

