

## Vaccine Carrier

Model No :- **GKS-VC-43**

**Tested & Approved by WHO accredited lab as per WHO/PQS/E004/VC01-VP.2**

### Specifications

Vaccine Storage capacity As per PQS	1.35 Liters
Gross Storage capacity	1.35 Liters
Weight Fully Loaded	4.36 Kgs
Weight- Empty (with empty ice packs)	2.22 kgs.
External material	HDPE
Internal lining material	HIPS
Insulation material	CFC Free Polyurethane
Insulation thickness	35 to 40 mm
External dimensions(LxWxH)	24x24x30 cm
Vaccine storage dimensions	9.9x9.9x16.5 cm
Lid type & Fixing	Removable
Number of ice packs supplied	4
Ice-Pack type	0.3 l compliant with PQS/E05/IP1.1
Volume of Ice-Pack	0.3 Liters
Cold Life without opening at +43°C	33 Hours
Warm life at -23°C	15 Hours
Cool life +43°C	7 Hours



### Address

**GKS Healthsol LLP**  
 B-1/03, Palm Grove Villas,  
 Ardee City, Sector 52, Gurugram,  
 Haryana 122011

### Factory

Plot 362, Sector 8 IMT Manesar,  
 Gurugram, Haryana 122050

### Contact us

email: [dinesh@gkshealth.com](mailto:dinesh@gkshealth.com)  
 Mob.No. +91 9899632558

Visit us:  
[www.gkshealth.com](http://www.gkshealth.com)