

Preserving
Vaccines,
Protecting
Lives



Vaccine Carrier

MODEL: GKS-VC-43

Storage
Capacity
1.44 Litres

Weight
Fully-Loaded
4.3 Kgs

Cold Life
35 Hrs &
46 Mins

Warm Life
16 Hrs &
56 Mins

Vaccine Carrier

MODEL: GKS-VC-43



PRODUCT SPECIFICATIONS

Storage Capacity	1.44 Litres	Warm Life	16 Hrs 56 Mins
Weight Fully-Loaded	4.3 Kgs	Cold Life At +43°C	35 Hrs 46 Mins
External Dimensions	24.5 x 24.5 x 29 cm	No. Of Ice Packs	4
Insulation Thickness	38 mm	Ice-Pack Type	0.3 Litres
WHO PQS Reference	WHO/PQS/E004/VC01-VP.2		