

Preserving
Vaccines,
Protecting
Lives



Vaccine Carrier

MODEL: GKS-VC-44

Storage
Capacity
1.39 Litres

Weight
Fully-Loaded
4.5 Kgs

Cold Life
41 Hrs

Warm Life
16 Hrs &
28 Mins

Vaccine Carrier

MODEL: GKS-VC-44



PRODUCT SPECIFICATIONS

Storage Capacity	1.39 Litres	Warm Life	16 Hrs 28 Mins
Weight Fully-Loaded	4.5 Kgs	Cold Life At +43°C	41 Hrs
External Dimensions	24.5 x 24.5 x 29 cm	No. Of Ice Packs	4
Insulation Thickness	38 mm	Ice-Pack Type	0.4 Litres
WHO PQS Reference	WHO/PQS/E004/VC01-VP.2		