

Preserving  
Vaccines,  
Protecting  
Lives



# Freeze-Free Vaccine Carrier

MODEL: GKS-FFVC-44-LR Grade A

Storage  
Capacity  
1.7 Litres

Weight  
Fully-Loaded  
6.4 Kgs

Cold Life  
30 Hrs &  
22 Mins

Cool  
Down Time  
64 Mins

# Freeze-Free Vaccine Carrier

MODEL: GKS-FFVC-44-LR



## PRODUCT SPECIFICATIONS

Storage Capacity	1.7 Litres	Cool Down Time	64 Mins
Weight Fully-Loaded	6.4 Kgs	Cold Life At +43°C	30 Hrs 22 Mins Long Range
External Dimensions	28.5 x 30.5 x 29.5 cm	No. Of Ice Packs	4
Insulation Thickness	40 mm	Ice-Pack Type	0.4 Litres
WHO PQS Code	E004 / 071	WHO PQS Reference	WHO/PQS/E004/VC02.1