

Preserving
Vaccines,
Protecting
Lives



Freeze-Free Vaccine Carrier

MODEL: GKS-FFVC-44-SR Grade A

Storage
Capacity
1.17 Litres

Weight
Fully-Loaded
5.9 Kgs

Cold Life
27 Hrs &
02 Mins

Freeze-Free Vaccine Carrier

MODEL: GKS-FFVC-44-SR



PRODUCT SPECIFICATIONS

Storage Capacity	1.17 Litres	Cold Life At +43°C	27 Hrs 02 Mins Short Range
Weight Fully-Loaded	5.9 Kgs	No. Of Ice Packs	4
External Dimensions	29 x 28.5 x 278 cm	Ice-Pack Type	0.4 Litres
Insulation Thickness	32 mm		
WHO PQS Reference	WHO/PQS/E004/VC02.1		